



www.healthcare.philips.com

Final Site Preparation Support Document

The equipment components shown in this drawing package are based on the current proposed purchase and are subject to change if modifications are made to the configuration.

VA Medical Center
Gainesville
Gainesville, FL



*Photo shown is not site specific.

Table of Contents	
Sheet Number	Sheet Name
CS	COVER SHEET
GN	GENERAL NOTES
A1	ARCHITECTURAL PLANS
AD1	EQUIPMENT DETAILS
AD2	EQUIPMENT DETAILS
S1	STRUCTURAL PLANS
S2	STRUCTURAL PLANS
SD	STRUCTURAL DETAILS
E1	ELECTRICAL PLANS
ED	ELECTRICAL DETAILS
NN	NETWORK NOTES
R	REFERENCE 11x17

Important Note:					
THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS. Philips assumes no liability nor offers any warranty for the fitness or adequacy of the premises or the utilities available at the premises in which the equipment is to be installed, used, or stored.					
For Architects and/or Contractors: The latest revision listed must be thoroughly reviewed so that all changes can be incorporated into your project, final revisions are valid for 90 days from latest revision date. Please check with Philips PM for current updates.					
Rev.	Date	Revision Descriptions	Planner	CPM	Approved By
-	8/18/2025	Created Preliminary Site Support Document per Quote #: Q-00552854. System upgraded to Rel. 3 from Rel. 2.2 reference project #: N-SOU210414. Created drawing packet in reference to project #: N-SOU070960 for existing conditions.	Pablo Arosemena	Michelle Horvath	
A	8/26/2025	Created Pre-Order Site Support Document per Quote #: Q-00552854.	Nicole Romero	Michelle Horvath	
B	11/14/2025	Created Final Site Preparation Support Document per Order #: 6600762648.	Andrew Defelice	Michelle Horvath	
C	12/4/2025	Replaced the Styker Booms with Drager Booms per CAD file: F02 VA Malcom Randall EP Lab_11NOV2025 AP.dwg-matched Styker locations per request. Updated Intrasight to go to EP Boom. Removed TV2 & TV3. Added Xper Hemo. Relocated AVOX2 & AVOX3. Updated conduits and added WR4. Moved EPB 0'-8" plan south. Reconfigured equipment room including IC and SH. Updated Planning Issue Notes.	Alison Baillargeon	Alan Rooks	

NOTE:
ALL SHEETS IN THE "R" SERIES
ARE REFERENCES TO VIEWS IN
THE MAIN SET AND ARE MEANT
TO BE PLOTTED AT 11x17

VA Medical Center
Gainesville
Gainesville, FL
PROVISIONAL

Drawing Title
COVER SHEET
THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.

Project
Azurion 7 B20/12 (Cardiac) - Swivel Catalyst
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com
Drawn By: Alison Baillargeon

Project Details
Drawing Number
N-SOU250567_RevC
Date Drawn: 8/26/2025
Quote: Q-00552854
Order: 6600762648

CS

© Koninklijke Philips Electronics N.V. 2019. All rights reserved. Reproduction in whole or in part is prohibited without prior written consent of the copyright holder.

HVAC Requirement for General Equipment Locations	
These temperature and humidity levels must be maintained in all (3) rooms (equipment, examination and control rooms).	
Temperature	59°F (15°C) to 86°F (30°C)
Temperature gradient	Max. 1°F / Minute (0.5°C / Minute)
Humidity (non-condensing) Humidity shall be stable within 10%	20% to 80%
Exam Room	*6479 BTU/hr
Equipment Room	*13299 BTU/hr
Control Room	*1876 BTU/hr
*Average heat emission during clinical use Data applicable for basic system: Large monitor + 4 x small monitor in Monitor Ceiling Suspension 1 workstation + 2 x small monitor in Control Room Add 8750 BTU/hr for optional Soomec UPS Add 1075 BTU/hr for additional large monitor Add 256 BTU/hr for additional small monitor Add 631 BTU/hr for additional workstation See AD details for additional heat load in case of UPS. Equipment's designed airflow is from front/side to back. Please design the air handling in the rack cabinet equipment area accordingly. Philips IGT-S systems may only be operated at an altitude of max. 9,843 ft (3,000m) above sea level as defined in our instructions for use (IFU).	
(25.0)	
Electrical Requirements Mains MA Cabinet	
Maximum Rated Power: 100KW	
Supply Configuration:	3 phase, equally sized insulated power conductors and an insulated equipment grounding conductor. Insulated grounding conductor shall have the same or larger size than line conductors. Line wires shall be no smaller than 4 AWG, 90°C or higher temperature rating. The conductor size is dependant on the upstream circuit breaker rating.
Nominal Line Voltage:	380 (Canada) -480 VAC, 60 Hz
Branch Power Requirement:	100 kVA (System only; verify UPS power requirements)
System Circuit Breaker:	3 phase, Type D with long-time delay and shunt trip
System Circuit Breaker: 80A rating. UPS breaker: 125A-150A, pending configuration.	
(24.1)	
Remote Control of Room Lighting	
Room lighting controls are the responsibility of the customer. Refer to sheet 3-ED for X-ray in use and room light switching requirements.	
(25.0)	
Milestone and General Conditions for Successful IGT Projects	
1. Planning, Design, and Implementation Process The multidisciplinary project team should be assembled as early as possible in the design process. The multidisciplinary team should include administrators, clinicians, infection preventionists, architects and other design professionals, facility managers, safety officers, security managers, users of equipment, and support staff relevant to the areas affected by the project as well as those with knowledge of the organization's functional goal for the project. Inclusion of patient advocates/consumers, A/E consultants, and construction specialists should be considered. Additionally, adherence to FGI (Facility Guidelines Institute) and ADA (Americans with Disabilities Act) guidelines should be taken into account.	
2. Schedule The customer or general contractor shall provide Philips with a project/construction schedule with milestones to assist in the coordination of delivery of Philips supplied products and primary equipment. Project schedule must be provided by customer to obtain production slot and delivery.	
3. Responsibility The customer shall be solely responsible, at their expense for preparation of site. Philips required specifications and any required MEP, construction and structural alterations shall be incorporated into customer's design and construction documents. Compliance with all safety, electrical, and building design codes relevant to the build out of the clinical area for Philips equipment and its installation is the customer's responsibility. Sufficiency of such plans and specifications, specifically including, but not limited to the accuracy of the dimensions described therein, shall be the sole responsibility of the customer. The customer shall advise Philips of conditions at or near the site, which could adversely affect the function of the equipment and/or carrying out of the delivery and installation work.	
a. Customer's structural engineer shall provide Philips with written certification that structural supports meet Philips requirements to permit delivery and installation of equipment.	
b. Customer shall acknowledge the final site preparation confirmation document.	
c. Upon completion of project, Customer's Architect and Engineers of record shall provide a set of As-Built project construction documents (.dwg) to Philips for closure of the Philips project history file.	
This shall ensure that such conditions are compliant and that the site is fully prepared and available for Philips before the installation work is due to begin.	
4. Permits Customer shall obtain all permits and licenses required by federal, state/provincial or local authorities in connection with the construction, installation and operation of the products and shall bear any expense in obtaining same or in complying with any related rules, regulations, ordinances and statutes.	
5. Infection Control and Interim Life Safety Measure The customer shall provide all means and methods necessary for compliance with Infection Control (IC) and Interim Life Safety Measures (ILSM) in connection with the construction and installation/operation of the products shown herein and shall bear any expenses related to same.	
6. Radiation Protection The customer or their contractor, at their own expense, shall obtain the service of a licensed radiation physicist to specify radiation protection and testing.	
7. Asbestos and Other Toxic Substances Philips assumes that there is no hazardous material contained in the project site. The customer is responsible for the removal of any materials, including but not limited to asbestos, deemed hazardous by local authorities, the EPA, OSHA, or any other authority having jurisdiction over the work. If such materials are discovered at any time that the work is proceeding, the work will immediately cease, the owner will be notified, and the work will again proceed after the owner has removed all of the hazardous material from the job site.	
8. Labor In the event local labor conditions make it impossible or undesirable to use Philips' regular employees for such installation and connection, such work shall be performed by laborers supplied by the customer, or by an independent contractor chosen by the customer at the customer's expense, and in such case, Philips agrees to furnish adequate technical support and supervision for proper completion for the Philips equipment installation.	
9. Extended Installation or Turnkey Work by Philips If an extended installation or turnkey contract exists between Philips and the customer for room preparation some of the responsibilities of the customer as depicted in these drawings may be assumed by Philips. In the event of a conflict between the work described in the turnkey contract work scope and these drawings, the turnkey contract work scope shall govern.	
10. Boom Vendors/3rd party cables Separation of cables shall be provided by boom vendor for the installation of cables within the boom. For example, use liquid-tight plastic conduit.	
(25.0)	

Customer Project Manager - Site Readiness Checklist				(25.0)	
SITE PREP TO BE COMPLETED PRIOR TO FIRST MILESTONE DATE / REQ/ PRIOR TO DELIVERY					
Site Planning Requirements/Checklist		REQUIRED 15 DAYS PRIOR TO FINAL PHILIPS INSPECTION		Room: cabinets, furniture, millwork, and counter tops installed.	
System Orientation: Verified per Philips Final Drawings		Customer (contractor) completed the Site ready validation (Unistrut conformance) email acknowledgment. Document archived in note/attachments.		Network Connections: Hardware (J-boxes, conduits, and network cable) is installed per Philips Final Drawings.	
FlexArm Clearance (if applicable): Verified per Philips Final Drawings.		Clearances: Verified to the closest obstacles (i.e. walls, cabinets), in order to lift up the C-arm, monitor support, etc.		UPS: Fully installed per Philips Final Drawings, and startup has been scheduled with vendor.	
3rd Party Booms (if applicable): Compatibility, locations, and isolation kits have been verified with boom vendor and the Philips Final Drawings are updated with the information. If required, a 'Request for Modification' has been submitted.		Cable Trough/Cable Raceway/Conduit: Installed, cleaned and locations per Philips Final Drawings. Cable openings are clear, without sharp edges.		UPS: Comm cable is run from UPS to M Cabinet	
Video Connection Boxes: Locations, video sources, and display destinations are verified with customer, and the Philips Final Drawings are updated with the information.		Conduit worksheet documenting validation of electrician's run length measurements has been uploaded to PSA.		HVAC (Climate Equipment): Customer/Contractor has verified system is operational and balanced to meet humidity and temperature requirements per Philips Final Drawings.	
Med Gas Box (if applicable): Location does not interfere with the installation and movement of the Philips equipment. The Philips Final Drawings are updated with the location.		Greenlee measuring tape (part no. 435 or equal) installed and validates run lengths.		Door Interlock Switch: If required, is installed per Philips Final Drawings.	
Drawings (Final): Approved by Customer shows all room obstacles to include millwork, lighting overlay, structure overlay, medical gases and plumbing.		Mains Power Supply: Installed per Philips Final Drawings. (Including impedance,isolated grounds, wire size verified, and distribution unit has been installed).		Customer Site Preparation: Verified per Philips Final Drawings.	
SITE PREP TO BE COMPLETED PRIOR TO FIRST MILESTONE DATE / REQUIRED PRIOR TO DELIVER		Back Boxes: Check label that the back boxes are manufactured by The DMG Group, installed with required covers and grommet material per Philips Final Drawings. Specifically, the spacing between the boxes and height off of finished floor.		PRE-DELIVERY DAY REQUIREMENTS	
GC Schedule aligns with Philips equipment delivery date. Permits and Inspections: completed by applicable governing authorities		Ceiling Plate for Equipment Rack (EP Boom) (if applicable): Installed and leveled per Philips Final Drawings.		3rd Party Booms (if applicable): Installed prior to Philips equipment delivery.	
Performance Testing Requirements Identified: Determine if certificate of compliance is required		REQUIRED 10 DAYS PRIOR TO FINAL PHILIPS INSPECTION		Installation Team: Has received the room drawings and necessary contact phone numbers.	
Existing Equipment: Is dismantled and removed from the site		Hospital Mains Supply: Allura R8.2 and Azurion, Installed per Philips Final Drawings, for connection in Cabinet Rear Cover (CRC) of MA-Cabinet.		Philips Project Space: Is clean, free of dust, all construction-related debris and tools have been removed.	
REQUIRED 30 DAYS PRIOR TO FINAL PHILIPS INSPECTION		Hospital Mains Supply: For system versions prior to Allura R8.2 (in the USA), the hospital mains must be available for connection at the PDU or gssPDU, and then to MA & ME-Cabinets (by the electrician).		Site Is Safe To Work: PPE requirements identified (Construction and Hospital). No open mains, slippery floors, sharp edges, or hazardous goods on site.	
Ceiling Unistruts (P1001 or engineer specified): Installed and leveled per Philips Final Drawings.		ERB Grounding Block: Installed within 6 feet or 150cm of building steel, and non-Philips-provided electrical equipment ground conductors, to include wall outlets installed per Philips Final Drawings.		Delivery Path and Truck Parking: Has been checked with the customer, including verifying floor loading, delivery route, elevator capacity, height, width and depth clearances, and a plan for bad weather.	
Fixing Blocks: Provided by Philips, properly seat in the Unistrut channel with no obstructions, as designed.		ERB Grounding Block: All Philips-provided electrical boxes and contractor-provided raceway are bond grounded to the ERB.		Transport & Handling Tools: Crane, forklift, wheels and trolleys have been specified with the LMP/rigging company. NOTE: If rigging provided by Philips, verify the vendor is on the Philips' approved suppliers list.	
Unistrut Height: Verified per Philips Final Drawings. Single plane measured from finished floor to bottom of Unistrut. For Bi-Plane, from Clea plate to bottom of Unistrut. Measurements have been documented and archived within PSA "Notes and Attachments".		Ceiling (False): Installed.		Site Access: Is available for after hours. Storage for tools, parts, covers and packing material has been arranged.	
Unistrut is square to the edge of the floor plate(s).		Ceiling (Hard): Installed		Room Security: Room is secure, with keys and alarm codes provided	
Unistrut validation workbook uploaded to PSC		Room Lighting: Installed and operational.		Parking: Parking area identified for installers.	
Floor Plates: Patient support and stand (if applicable), are installed, isolated, and leveled, at the correct locations per Philips Final Drawings. For Biplane systems, the plates are spaced correctly and squared with each other.		Restricted ceiling area "A" is free from any objects that project below the finished ceiling.		Restroom Facilities: Toilet facilities, including area to wash up, are available.	
Floor Levelness: Checked with laser level and is level per Philips Final Drawings.		Restricted ceiling Area "B" is free of any objects that project more than the specified maximum distance below the finished ceiling.		CPM has coordinated any ongoing work that will overlap with the installation window with the lead installer	
Table Isocenter: Verified per Philips Final Drawings.		Ceiling Obstructions: Verify there are no obstructions where Philips rails will be installed.		Philips Field Pre-Job Safety Brief completed and posted on door.	
Walls: Installed.		REQUIRED 5 DAYS PRIOR TO FINAL PHILIPS INSPECTION		PRIOR TO PHILIPS SYSTEM POWER UP	
Sprinklers: Installed.		Lead-Lined Doors: Installed.		Wall Outlets: Installed and functional.	
SP and TV J-Boxes: Located per Philips Final Drawings.		Lead-Shielded Glass: Installed.		PRIOR TO INSTALL COMPLETE	
CPM has updated PSC with facility IT contact information		Walls: Final finished, (i.e. final coat painted and/or tiled).		Physicist: If required, verify the Physicist has been scheduled.	
Video Wall Boxes: Located per Philips Final Drawings		Hard Ceiling: Painted.		UPS: Commissioned and certified by UPS vendor.	
HVAC (Climate Equipment): Installed.		Flooring: Installed and covered with protective covering (i.e. scratch protection).		X-Ray in Use or Warning Light: If required, is installed per Philips Final Drawings.	
REQUIRED SITE READINESS - TABLE OF CONTENTS					(25.0)
Items to be completed Prior to Equipment Delivery:		Inspected by:	Date:	Common Site Drawing Reference: Sheet/Section/Detail	
Table Iso Center: Verified per Philips Final Drawings.				1-A1,2-A1.	
Floor Levelness: Checked with Laser Level.				1-S1, S1/Equipment Support Information.	
Floor Plates: Installed, isolated, and leveled at the correct locations.				1-S1, 2-S1, 3-S1, S1/Equipment Support Information.	
Ceiling Unistruts (P1001): Installed and leveled per Sheet with written engineer certification.				2-S2, 4-S2, 6-S2, S2/Equipment Support Information.	
Modular Laminar Array with Unistruts (P1001) installed and level per sheet if Class 3 required.				4-S2, S2/Equipment Support Information.	
Ceiling Plate for Equipment Rack (EP Boom) (if applicable): Installed and leveled.					
Cable Trough/Raceway/Conduit: A Greenlee measuring tape (part no. 435) are Installed, verified and locations checked.				E1/Point to Point Run Lengths, 1-E1,4-E1.	
3rd Party Booms (if applicable): Structural support installed and verified with boom vendor and locations.					
Ceiling Height: Verified and measure from bottom of Unistrut.				3-S2, 4-S2.	
Verify the Fixing Blocks sits properly in the Unistrut channel with no obstructions.				4-S2,1-S2.	
Ceiling Obstructions: Verify there are no obstructions where Philips rails will be installed.				1-S2, 3-S2.	
Clearances: Verified to the closest obstacles (i.e. walls, cabinets), in order to lift up the C-arm, monitor support, etc.				1-AD1.	
Back Boxes: Installed with required covers and grommet material.				1-E1, 4-ED, 5-ED.	
ERB Conductor Bar Installed per drawing.				1-E1,2-E1.	
All electrical boxes and raceway are grounded to the ERB.				1-ED, 2-ED, 5-ED	
Mains Power Supply: Installed per drawing (Including impedance, isolated grounds, wire size and circuit breakers verified).				ED/Power Quality Requirements (Azurion), ED/Branch Circuit and Wire Gauge Requirement (Azurion).	
Mains Supply Wiring installed for connection in Cabinet Rear Cover (CRC) of MA-Cabinet.				5-ED	
Video Connection Boxes: video sources, and display destinations are verified with customer and located.				A1/AUX Chart.	
Med Gas Box (if applicable): Location does not interfere with the installation and movement of table.				3-S1	
Walls: Installed and final finished per customers architectural specifications.					
Millwork: Completely installed in all rooms.					
Leaded glass installed.					
Flooring: Installed and covered with protective covering (1/8" masonite).				S1/Equipment Support Information.	
Lighting and Wall Outlets Installed and functional.				ED/General Electrical Information.	
HVAC system is operational and commissioned to engineered specifications.				GN/HVAC Requirement for General Equipment Locations.	
3rd Party Booms (if applicable): Installed prior to Philips equipment delivery.					
UPS: Fully installed per Philips Final Drawings, and startup has been scheduled with vendor.				5-ED.	
X-Ray in Use Light is installed.				3-ED.	
Physicist: If required has been scheduled.					
UPS: Commissioned and certified by UPS vendor.					
Permits Inspections: Completed or scheduled by applicable governing authorities.					
All network information provided by facility IT, i.e. IP addresses (static IPs only), AE Titles, SNM, GTWY and DNS server are available.				NN.	
Project Space: Is clean, free of dust, all construction-related debris and tools have been removed.					
PPE requirements identified (Construction and Hospital). No open electrical or hazardous materials on site.				GN/Milestone and General Conditions for Successful IGT project.	
Site Access Is available for after hours. Storage for tools, parts, covers and packing material has been arranged.					
Room is secure with keys provided.					
Complete Delivery Path Route and backup for bad weather has been reviewed & reverified with the customer and lead FSE.					

INSPECTIONS		
(24.0)		
	Inspected by:	Date:
Facilities:		
Clinical:		
Infectious Control:		
Architect:		
Engineer:		
Administrator:		
Project Manager:		
Sales Specialist:		
Account Manager:		
Service:		

AVO NAMING	
Line Item	Corresponding Naming
FCV0972 (powered by system, in control room)	AVOC
FCV0973 (powered by system, in exam room)	AVOE
FCV0974 (powered by hospital)	AVOX
FCV0975 (4K AVO)	AVOK
FCV0976 (XL)	AVOL
Integration Kits	AVOB
FlexVision Non-Philips Frame	AVOM


PHILIPS

VA Medical Center
Gainesville

Gainesville, FL

PROVISIONAL

Drawing Title

GENERAL NOTES

THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.

Project

Azurion 7 B2012 (Cardiac) - Swivel Catalyst

VA Medical Center Gainesville

Gainesville, FL

Room: C246-1

Philips Contacts

Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com

Drawn By: Alison Baillargeon

Project Details

N-SOU250567_RevC

Date Drawn: 8/26/2025


Quote: Q-00552854

Order: 6600762648

GN

5 TRANSPORT DETAILS - CLEA FLOOR


AD1 SCALE: N.T.S.



(23.0)

6 TRANSPORT DETAILS - LATERAL ARC N


AD1 SCALE: N.T.S.



(24.0)

7 MONITOR SUSPENSION DETAIL

AD1 SCALE: N.T.S.




AD1

SCALE: N.T.S.


1 SYSTEM DETAIL

AD1 SCALE: N.T.S.




2 SYSTEM CABINET DETAILS

AD1 SCALE: N.T.S.



4 DOCUMENTATION BOX DETAIL

AD1 SCALE: N.T.S.



Drawing Title

Project
Azurion 7 B20/12 (Cardiac) - Swivel Catalyst
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com

Drawn By: Alison Baillargeon


Project Details
Drawing Number
N-SOU250567_RevC
Date Drawn: 8/26/2025
Quote: Q-00552854
Order: 6600762648

EQUIPMENT DETAILS


VA Medical Center
Gainesville
Gainesville, FL
PROVISIONAL




THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.



1 SERVICE HUB
AD2 SCALE: N.T.S.




2 CONTROL ROOM EQUIPMENT ADAPTER
AD2 SCALE: N.T.S.




3 AUXILIARY EQUIPMENT ADAPTER
AD2 SCALE: N.T.S.

Applicable for all types of AVOX (AVON / AVOX / AVOK / AVOL / AVOB / AVOM)




4 AUXILIARY VIDEO OUTPUT DETAIL
AD2 SCALE: N.T.S.



Required Items	Configuration
Video Connection Box mounting and Video Connection Box-TX cabling	Present
Flexible conduit with guide wire	Present
Number of shelves	5
Number of drawers	1
Cable management accessories: 2 rails, 8 hooks	Present
Number of infusion bottle holder (1-pole)	1
Number of extension arms with VESA 100 monitor interface	1
Number of power outlets	4 standard and 4 emergency
Number of Ethernet network outlets	1

Note: This does not include the Video Connection Box itself, only a mechanical solution to mount a Video Connection Box. For details, consult the Sales Department.

*The above is standard configuration - not site specific.




- MED-Gas**
- The EP Boom comes with option MED-Gas outlets for customer convenience. Connection and interfacing of the MED-Gas to the EP Power Boom is the contractor's responsibility.
 - Standard configuration listed below. Other configurations may be possible. Contact Philips Project Manager for information.
 - 4 is standard and 8 is maximum in any configuration: Vacuum, N2O, Oxygen (or NO gases), WAGD
 - Chemtron, DISS, Puritan Bennett, or Ohmeda fitting type
 - 4 standard power duplex outlets
 - 4 emergency power duplex outlets
 - MAX total of 12 standard or emergency power duplexes of any combination
 - 2 network ports (2 minimum and 8 maximum)
 - 1 conduit in rack for 3rd party cabling (with 4 optional 1xG23 (22mm) and up to 4 max but not recommended to have more than 3 since went up to G23 from G17)
 - Equipment rack configuration must be confirmed by customer with Philips Project Manager prior to system entering production.
 - Once production starts, it is not possible to change the configuration of the equipment rack.



EPB	Equipment Rack DVI (EP Power Boom)	Weight
		865 lbs (393 kg)

Weight includes:
-Pre-installation ceiling frame 309 lbs (141 kg)
-Equipment Rack 389 lbs (177 kg)
-Accessories max. 167 lbs (76 kg)

5 EQUIPMENT RACK DVI (EP BOOM) DETAIL
AD2 SCALE: N.T.S.




Detail - Typical Control Room Layout with FlexSpot
(Not to scale - Not site specific)




7 TABLE MOUNTED INJECTOR, REMOTE PANEL AND CONSOLE
AD2 SCALE: N.T.S.

6 CONTROL ROOM JUNCTION BOX DETAIL
AD2 SCALE: N.T.S.

CY	Control Room Junction Box (All Components)	Weight
		86 lbs




8 INTRASIGHT - WORKSTATION DETAIL
AD2 SCALE: N.T.S.




UPS Cabinet	Socomec - Low Load Fluoro (LLF) UPS	Weight	Heat Dissipation
UPS		1873 lbs	8750 BTU/hr
SBO		6.7 lbs	- BTU/hr

Notes:
* Required ventilation for the front - 7 1/8" (200mm)
* Required ventilation for the rear - 7 1/8" (200mm)
* Required ventilation for the top - 6" (155mm)
* Acoustic noise level: <= 58 dB(A) @ 1 meter in front of the rack and 1 meter high (1 meter = 39.37")

9 SOCOMECE LLF UPS DETAIL
AD2 SCALE: N.T.S.



10 MAVIG 2.5m CEILING TRACK DETAIL
AD2 SCALE: N.T.S.



XHA	Xper Hemo App	Weight
		84 lbs
tUPS	Tripp Lite UPS	Weight
tUPS2		96.2 lbs
IX3	Intellivue X3	Weight
		- lbs

11 XPER HEMO APP
AD2 SCALE: N.T.S.



VA Medical Center
Gainesville
Gainesville, FL
PROVISIONAL

EQUIPMENT DETAILS
THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSIDERED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.

Project
Azurion 7 B2012 (Cardiac) - Swivel
Catalyst
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com

Drawing Title

Project Details
Drawing Number
N-SOU250567_RevC
Date Drawn: 8/26/2025
Quote: Q-00552854
Order: 6600762648

AD2

Equipment Support Information

1. General

The customer shall be solely responsible, at its expense, for preparation of the site, including any required structural alterations. The site preparation shall be in accordance with this plan and specifications, the architectural/construction drawings and in compliance with all safety and building codes. The customer shall be solely responsible for obtaining all construction permits from jurisdictional authority.

Customer's structural engineer shall provide Philips with written certification that structural supports meet Philips requirements to permit delivery and installation of equipment. Upon completion of project, Customer's Architect and Engineers of record shall provide a set of As-Built project construction documents (dwg) to Philips for closure of the Philips project history file.

2. Equipment Anchorage

Philips provides, with this plan and specifications, information relative to equipment size, weight, shape, anchoring hole locations and forces which may be exerted on anchoring fasteners. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings, information regarding the approved method of equipment anchoring to floors, wall and/or ceiling of the building. Any anchorage test required by local authority shall be the customer's responsibility. Stud type anchor bolts should not be specified as they hinder equipment removal for service. Consult with Philips service prior to specifying anchor methods. Philips equipment must be electrically isolated from anchorage.

3. Floor Loading and Surface

Philips provides, with this plan and specifications, information relative to size, weight and shape of floor mounted equipment. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings confirmation of the structural adequacy of the floor upon which the equipment will be placed. Any load test required by local authority, shall be the customer's responsibility.

The floor surface/base plates upon which Philips equipment is to be located shall be flat and level to within 0.05 degrees.

4. Seismic Anchorage (For Seismic Zones Only)

All seismic anchorage hardware, including brackets, backing plates, bolts, etc., shall be supplied and installed by the customer/contractor unless otherwise specified within the support legend on this sheet. Installation of electronic cabinets to meet seismic anchorage requirements must be accomplished using flush mounted expansion type anchor/bolt systems to facilitate the removal of a cabinet for maintenance. Do not use threaded rod/adhesive anchor systems. Philips equipment must be electrically isolated from anchorage. Seismic anchoring information are available upon request. Architect of Record to determine seismic planning needed for site. Consult with Philips regarding any anchor system issues.

5. Floor Obstructions/ Floor Coverings

There shall be no obstructions on the floor (sliding door tracks, etc.) within the serviceability area of the Philips technical cabinets. Floor must be clear to allow cabinets to be pulled away from the wall for service. Technical equipment room floor shall be commercial grade "VCT" Vinyl Composition Tile or a flooring material of equal hardness and compression resistance.

(25.0)

Important Notes:

1. Anchors for items that are installed/anchored by customer/contractor shall be provided by customer/contractor.
2. Anchors for items that are installed/anchored by Philips shall be provided by Philips. If customer's engineering documents specify anchors other than those listed in this document, the anchors shall be provided by customer/contractor and installed by Philips.
3. In all instances, the wall and/or floor support are the sole responsibility of the customer/contractor. The customer's architect/engineer of record shall specify wall and/or floor support sufficient for the bolt forces shown on the details.

Floor & Wall Support Legend

Floor & Wall Support Legend		
A. Furnished and installed by Philips B. Furnished by customer/contractor and installed by customer/contractor C. Installed by customer/contractor D. Furnished by Philips and installed by contractor E. Existing F. Future G. Optional H. Furnished by Philips and installed by Third Party		
Item Number	Description	Detail Sheet
D AD7	AD7 Swivel Floor Plate	S1/3
D CF	Floor Clea Floor Plate	S1/5
B SH	Service Hub	SD/3

Important Notes:

- * The Iso Center is the initial and most crucial control point in the construction of a procedure room. Once established, it guides the placement of all imaging equipment, third-party devices, mechanical, electrical, and plumbing systems (MEPs), as well as the overall construction. These are positioned relative to the Iso Center's X, Y, and Z coordinates to avoid operational issues and clinical limitations.

Equipment Support Information


11. Requirements for External Vibration

The maximum allowed external vibration level of floors and ceilings, to which the equipment is mounted that will not adversely affect the image quality, is specified in terms of RMS velocity levels in 1/3-octave or terts bands, as follows:

Center frequency [Hz]	1	1.25	1.6	2	2.5	3.15	4	5	6.3	8
Terts band value [$\mu\text{m/s}$] (RMS)	50.8	50.8	50.8	50.8	63.5	80.01	101.6	127	160	203.2
Center frequency [Hz]	10	12.5	16	20	25	31.5	40	50	63	80
Terts band value [$\mu\text{m/s}$] (RMS)	254	317.5	406.4	406.4	406.4	406.4	406.4	406.4	406.4	406.4
Center frequency [Hz]	100	125	160	200						
Terts band value [$\mu\text{m/s}$] (RMS)	406.4	406.4	406.4	406.4						

Terts Band Specification for External Vibrations

A graphical representation of this specification is given below:




Terts Band Specification for External Vibrations

Terts band spectra shall be calculated on the basis of time traces with a duration of 10 minutes (600 seconds), taken at representative locations and during representative times during working days.

(23.1)

Detail - Floor Structural Azurion 7 B20/12, B20/15 - Swivel



Floor plates supplied by Philips / installed by contractor.

Clea Floor	AD7 Swivel Table
Floor Plate to Floor Bolt Forces:	Floor Plate to Floor Bolt Forces:
(Tension) T _{max} = 2668 lbs/bolt	(Tension) T _{max} = 2332 lbs/bolt
(Shear) V _{max} = 892 lbs/bolt	(Shear) V _{max} = 1164 lbs/bolt

Azurion 7 B20/12, B20/15 - Swivel Structural Detail

2
S1

SCALE: N.T.S.


(24.0)

1
S1

FLOOR & WALL SUPPORT LAYOUT

SCALE: 1/4" = 1'-0"

Required Unistrut Height: 9' - 9 3/8" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.



VA Medical Center
Gainesville

Gainesville, FL

PROVISIONAL

Drawing Title

STRUCTURAL PLANS

THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.

Project
Azurion 7 B20/12 (Cardiac) - Swivel

Catalyst

VA Medical Center Gainesville

Gainesville, FL

Room: C246-1

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com


Drawn By: Alison Ballageon

Project Details
Drawing Number
N-SOU250567_RevC
Date Drawn: 8/26/2025
Quote: Q-00052854
Order: 6600762648

S1

Detail - AD7 SyncraTilt/ Tilt Table, Swivel Base - Clearance Area for Table in Normal Position

(Not to scale)



(23.0)

AD7 SWIVEL PLATE - TABLE CLEARANCE DETAIL

3
S1

SCALE: N.T.S.

Detail - AD7 Swivel Floor Plate - Notes for Installation



(Not to scale)

1. 1.18" (30mm) thick floor plate, flush mounted with top of slab.
2. Level within 1/32" (1mm) per 39" (1m) across surface of plate.
3. Anchor bolts with washer must not protrude above surface of floor plate.
4. Align center of AD7 swivel floor plate with center of Clea floor plate.

Floor plate mounting to the building:

In case threaded rods are used and the nut protrudes above the floor plate surface, DO NOT GRIND DOWN THE NUT, but follow the procedure stated below to ensure the nut is flush with the floor plate surface.

1. Use Jam nuts M16 (h=8.0mm) or (h=9.5mm).
2. Use only 1 washer.
3. Use locite 243 instead of a lock washer.
4. Torque value to be determined by structural engineer of record.





(25.0)

Detail - Floor Clea Floor Plate - Notes for Installation

(Not to scale)

1. 1 3/8" (20mm) thick floor plate, flush mounted with top of slab.
2. Level within 1/32" (1mm) per 39" (1m) across surface of plate.
3. Anchor bolts with washer must not protrude above surface of floor plate.
4. Clea floor plate reference line (Installation alignment holes) are 1/8" (10mm) offset from center of plate.
5. Do NOT use the center of Clea floor plate to align the floor plate with AD7 floor plate.
6. Align Clea floor plate reference line (Installation alignment holes from Clea floor plate; 315mm from edge of Clea floor plate), with the center of AD7 floor plate.




Note: Use all ten holes for anchoring the floor plate.
Ten holes, 0.75" (19mm) dia.
Counterbored 1.38" (35mm) dia. x 0.56" (14mm) deep
1'-2"

Floor plate mounting to the building:

In case threaded rods are used and the nut protrudes above the floor plate surface, DO NOT GRIND DOWN THE NUT, but follow the procedure stated below to ensure the nut is flush with the floor plate surface.

1. Use Jam nuts M16 (h=8.0mm) or (h=9.5mm).
2. Use only 1 washer.
3. Use locite 243 instead of a lock washer.
4. Torque value to be determined by structural engineer of record.



Floor plate mounting to the building:

In case threaded rods are used and the nut protrudes above the floor plate surface, DO NOT GRIND DOWN THE NUT, but follow the procedure stated below to ensure the nut is flush with the floor plate surface.

1. Use Jam nuts M16 (h=8.0mm) or (h=9.5mm).
2. Use only 1 washer.
3. Use locite 243 instead of a lock washer.
4. Torque value to be determined by structural engineer of record.

(25.1)

CLEA FLOOR - INSTALLATION DETAIL

5
S1

SCALE: N.T.S.

© Koninklijke Philips Electronics N.V. 2019. All rights reserved. Reproduction in whole or in part is prohibited without prior written consent of the copyright holder.

Equipment Support Information

1. General
The customer shall be solely responsible, at its expense, for preparation of the site, including any required structural alterations. The site preparation shall be in accordance with this plan and specifications, the architectural/construction drawings and in compliance with all safety and building codes. The customer shall be solely responsible for obtaining all construction permits from jurisdictional authority.

Customer's structural engineer shall provide Philips with written certification that structural supports meet Philips requirements to permit delivery and installation of equipment. Upon completion of project, Customer's Architect and Engineers of record shall provide a set of As-Built project construction documents (.dwg) to Philips for closure of the Philips project history file.

2. Equipment Anchorage
Philips provides, with this plan and specifications, information relative to equipment size, weight, shape, anchoring hole locations and forces which may be exerted on anchoring fasteners. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings, information regarding the approved method of equipment anchoring to floors, wall and/or ceiling of the building. Any anchorage test required by local authority shall be the customer's responsibility. Stud type anchor bolts should not be specified as they hinder equipment removal for service. Consult with Philips service prior to specifying anchor methods. Philips equipment must be electrically isolated from anchorage.

3. Ceiling Support Apparatus
a. Philips provides, with this plan and specifications, information relative to size, weight and shape of ceiling supported equipment. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings, information regarding the approved method of structural support apparatus, fasteners and anchorage to which Philips will attach equipment. Any anchorage and/or load test required by local authority shall be the customer's responsibility. Philips equipment must be electrically isolated from anchorage.
b. Contractor to clearly mark Philips equipment longitudinal centerline on bottom of each Unistrut support.
c. The structural Unistrut surface to which Philips equipment is to be attached, shall have horizontal equipment attachment surfaces parallel, square and level to within .230" (6mm) per entire span.
d. Any drilling and/or tapping of holes required to attach Philips equipment to the structural support apparatus shall be the responsibility of the customer.

4. Lighting
Luminaires shall be placed in such a position that they are not obscured by equipment or its movement, nor shall they interfere with Philips ceiling rails and equipment movement or otherwise adversely affect the equipment. Such luminaire locations shall be the sole responsibility of the customer.

5. Ceiling Obstructions
There shall be no obstructions that project below the finished ceiling in the area covered by ceiling suspended equipment travel. See detail 1/S2.

6. Safety Factors
Ceiling loads as mentioned in the PRD are worst case loads and excluding safety factors. Proper safety factors need to be applied by Design Professionals/Engineer of Record.

(24.0)

Detail - Philips Fixing Block for Philips Ceiling Rails (Clip Rails)


(Not to scale)

General Requirements:
1. Philips does not specify the overhead equipment support structure. Unistrut may or may not be used. If Unistrut is used, it is up to Unistrut and the structural engineer for the project to determine which of its products are appropriate for each project.
2. P1001 Unistrut is specified for the ceiling rail connection point.

Finished Ceiling Requirements:
1. Finished ceiling height shall be 1/2" (6mm) to 1/2" (13mm) above bottom of Unistrut.

Fixing Block Installation Requirements:
1. Nothing shall be attached to the Unistrut with any fastener that protrudes into the Unistrut which would interfere with positioning of the fixing block.
2. Fixing blocks for Philips ceiling rails (Clip rails) are designed to be installed in P1001Unistrut.
3. The inside of the Unistrut must be clear of obstructions.

Unistrut Requirements:
1. Welding Unistrut may warp Unistrut and deteriorate the structural integrity of the Unistrut. Consult the Structural Engineer of Record prior to welding any Unistrut.
2. If third-party structural laminar modular array vendors are to be used, please contact site planning.
3. For third-party structural laminar modular array vendors, the minimum rail fixing block Unistrut length must be 7-5/8" (194mm).




(25.1)


4 CLIP RAIL CROSS SECTION DETAIL

S2 SCALE: N.T.S.


Ceiling construction is responsibility of local structural engineer



Double nut supplied by Philips




Vertical Stiffness




Side View X


Ceiling construction is responsibility of local structural engineer



4" (101 mm)



Horizontal Stiffness



Side View Y

Stiffness Requirements of Ceiling:
1. Horizontal Stiffness: preferred 10,000,000 Newton/meter - 57.1 kbl/in, minimal 6,000,000 Newton/meter - 34.2 kbl/in
2. Vertical Stiffness: preferred 10,000,000 Newton/meter - 57.1 kbl/in, minimal 6,000,000 Newton/meter - 34.2 kbl/in
3. Rotation Stiffness: minimal 100,000 Newtonmeter/Rad - 885 (kbl in)/Rad
4. Ceiling interface stiffness requirements at each ceiling connection point

For Clea Stand (1160 kg or 11600 N): the maximum allowed deflection of the ceiling and/or floor construction is 1.93 mm.
For Larc Stand (835 kg or 8350 N): the maximum allowed deflection of the ceiling construction is 1.39 mm.

(25.0)

5 CLIP RAIL STIFFNESS DETAIL

S2 SCALE: N.T.S.

Ceiling Support Legend

A Furnished and installed by Philips
B Furnished by customer/contractor and installed by customer/contractor
C Installed by customer/contractor
D Furnished by Philips and installed by contractor
E Existing
F Future
G Optional
H Furnished by Philips and installed by Third Party

Item Number	Description	
LR	2 - Philips Larc C Cardiac Rails	S2
MSR	2 - Philips Monitor Equipment Rails	S2
UNI	Unistrut (P1000/P1001 in meeting Philips ceiling requirements, geometry of channel and geometry of fixing block) - Bottom of Unistrut 1/4" (6mm) to 1/2" (13mm) Below Finished Ceiling	S2
MAV	Mavig 2.5m Ceiling Track 360	AD2
EPB	Equipment Rack DVI (EP Power Boom) Ceiling Plate and Interface Plate	SD
AO	Additional Opening on Finished Ceiling - Required for Installation of the Boom (to be covered after Equipment Rack DVI (EP Power Boom) is installed)	SD


Important Notes:

- The Iso Center is the initial and most crucial control point in the construction of a procedure room. Once established, it guides the placement of all imaging equipment, third-party devices, mechanical, electrical, and plumbing systems (MEPs), as well as the overall construction. These are positioned relative to the Iso Center's X, Y, and Z coordinates to avoid operational issues and clinical limitations.
- All Unistruts are to be true and square to the edge of the floor plate.

Planning Issues and Considerations

- Structural Engineer to verify reusability of Unistrut on site. The permissible stress and strain amount on existing Unistrut may be insufficient.
- Monitor rails were moved 2-3/8" plan west to meet minimum overhang requirements.
- System rails exceed minimum overhang (5-7/8").
- Architect to modify reflected ceiling plan to accommodate system ceiling equipment, Unistrut, and electrical ceiling boxes.
- Boom locations match Stryker locations. Draeger to update plans. Final boom placement to be discussed between customer and 3rd Party Vendor. Boom B1 moved 8" plan east to avoid collision with column and Boom A1 moved 1" plan north per request.

Detail - Restricted Ceiling Area for Objects that Project Below Finished Ceiling
(Not site specific)




(24.0)

3 RESTRICTED AREA DETAIL

S2 SCALE: N.T.S.


Detail - Ceiling Structural Azurion 7 B20/12, B20/15 - AD7 Swivel



(25.0)

AZURION 7 B20/12, B20/15 - SWIVEL STRUCTURAL DETAIL


S1 SCALE: N.T.S.



2 CEILING SUPPORT LAYOUT - UNISTRUT

S2 SCALE: 1/4" = 1'-0"

Required Unistrut Height: 9' - 9 3/8" + 3/8" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.



1 CEILING SUPPORT LAYOUT - EQUIPMENT

S2 SCALE: 1/4" = 1'-0"

Required Unistrut Height: 9' - 9 3/8" + 1/4" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.



VA Medical Center
Gainesville
Gainesville, FL
PROVISIONAL

STRUCTURAL PLANS
THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.

Project
Azurion 7 B20/12 (Cardiac) - Swivel
Catalyst
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com

Drawn By: Alison Ballargeon

Project Details
Drawing Number
N-SOU250567_RevC
Date Drawn: 8/26/2025
Quote: Q-00552854
Order: 6600762648

S2


Pre-Evaluated and -Approved Anchor Reference List for Philips Installers

Anchors for items that are installed/anchored by customer/contractor shall be provided by customer/contractor. Anchors for items that are installed/anchored by Philips shall be provided by Philips. If customer's engineering documents specify anchors other than those listed below, the anchors shall be provided by customer/contractor and installed by Philips. In all instances, the wall and/or floor support are the sole responsibility of the customer/contractor. The customer's architect/engineer of record shall specify wall and/or floor support sufficient for the bolt forces shown on the details.

Equipment	Option	Anchor Style (provided by Philips)	Anchor Size (provided by Philips)	Qty.	Support Size & Material (provided & installed by customer/contractor)
Mavig Ceiling Track	A	Bolts, flat washer, lock washer, spring nuts	A307 Grade or ASME Grade 5 Bolts: $\frac{3}{8}$ " (10mm) x 2" (50mm) L Spring Nuts: $\frac{3}{8}$ " (10mm)	8	Unistrut


(12.0)

Detail - Example Control Room Junction Box - Cable Outlets
(Not to scale)




(23.0)

Detail - Junction Box Support Frame Option
(Not to scale)



(23.0)


Detail - Control Room Junction Box - Hole Pattern For Mounting
(Not to scale)



(25.0)

1 CONTROL ROOM JUNCTION BOX MOUNTING DETAILS & ANCHOR REFERENCE LIST
SD SCALE: N.T.S.


Detail - AUXEA/CREA/AVO - Hole Pattern For Mounting
(Not to scale)



Applicable for all types of AVOs (AVOC / AVOE / AVOK / AVOK / AVOL / AVOB / AVOM)

(24.1)


2 AUXEA/CREA/AVO MOUNTING DETAIL
SD SCALE: N.T.S.



(23.0)

3 SERVICE HUB MOUNTING DETAIL
SD SCALE: N.T.S.


Detail - Calculation of Bolt Length Interface Plate Power Boom
(Not to scale)



(25.1)


Detail - Exploded Ceiling Substructure for Power Boom
(Not to scale)

To be installed by customer/contractor per engineer of record.
Substructure must be ordered separately. Contact Ondal Customer Service, customerservice@ondal.com.




(25.0)

Detail - EP Power Boom Round Interface Plate
(Not to scale)



Max. vertical force = 1060 lbs (481 kg)(4718N) (21.0)

Detail - Exploded Interface Plate Mount Power Boom
(Not to scale)



(21.0)

4 EQUIPMENT RACK DVI - ROUND PLATE DETAILS
SD SCALE: N.T.S.



VA Medical Center
Gainesville
Gainesville, FL
PROVISIONAL

STRUCTURAL DETAILS

THE INFORMATION IN THIS PACKAGE IS PROVIDED AS PHILIPS EQUIPMENT REQUIREMENTS, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS.

Drawing Title

Project
Azurion 7 B20/12 (Cardiac) - Swivel Catalyst
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

Philips Contacts

Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com

Drawn By: Alison Baillargeon


Project Details

Drawing Number
N-SOU250567_RevC
Date Drawn: 8/26/2025
Quote: Q-00552854
Order: 6600762648

SD

Definition & Abbreviations	
Terms	Description
*	"Mandatory items to be filled in when available in the Cathlab. Missing these, might impact the installation time."
AE title	An AE Title is used by an Application Entity (AE) to identify itself. AE Titles need to be locally unique and are typically managed by a system administrator.
BioMed	BioMedical Engineer
CMVPN	Customer Managed Virtual Practice
CT	Computed Tomography
DHCP	Dynamic Host Configuration Protocol
DICOM	Digital Imaging and Communication in Medicine
DNS	Domain Name System
FORS	Fiber Optic RealShape
FQDN	A fully qualified domain name (FQDN) is the complete domain name for a specific computer, or host, on the internet. The FQDN consists of two parts: the hostname and the domain name. For some registration part, this methodology can be used for registration.
FSE	Field Service Engineer
HCF	HealthCare Facility (= customer)
Hostname	In computer networking, a hostname is a label that is assigned to a device connected to a computer network and that is used to identify the device in various forms of electronic communication, such as the World Wide Web.
IP address	Packaging configuration containing one or more serviceable items. Optionally containing tools to install the serviceable items.
IWS	Interventional Workspot
MAC Address	A media access control address (MAC address) is a unique identifier assigned to a network interface controller (NIC) for use as a network address in communications within a network segment. MAC address can mostly be found on the rear side of the PC, or can be found by the survey in the service application
MPPS	Modality Performed Procedure Step
MR	Magnetic Resonance
NCS	Network Connectivity Sheet
NTP	Network Time Protocol
PACS	Picture Archiving and Communication System
PM	Project Manager, responsible for managing the cathlab and X-ray installation.
Port number	A port number is a way to identify a specific process to which an Internet or other network message is to be forwarded when it arrives at a server.
PRS	Philips Remote Services
RIS	Radiology Information System
RSE	Remote Service Engineer
RSg	Remote Services Gateway
RSGoP (Variant 1)	Remote Services Gateway on Premises using a Philips HS Edge
RSGoP (Variant2)	Remote Services Gateway on Cloud using a CMVPN
RSGoP (Variant 3)	Remote Services Gateway on Cloud using Intemt
SRS	Structured Report Server
WLM	Work List Management

IMPORTANT NOTE: It is the customer's responsibility to coordinate with the local Philips Engineer to provide ALL required network information and install ALL required network cabling & drops according to Philips specifications PRIOR to the scheduled installation start date. Failure to do so may delay system installation and jeopardize the customer hand over date.



Hospital Information			
Health Care Facility (HCF)		HCF - CathLab identification X-ray system	Responsibility
Hospital name	E.g. Hospital X	HCF IT	Cathlab / room name
Country		HCF IT	Cathlab/Room ID
City		HCF IT	Physical location in HCF
Postal code		HCF IT	
Address		HCF IT	
Philips Field Service Engineer		HCF Project Manager/ BioMed	Responsibility
Contact name	Philips PM	Email address	Philips PM
Phone number	Philips PM	Secondary contact name	Philips PM
Email address	Philips PM	Secondary contact Email	Philips PM
Phone number	Philips PM	Secondary contact Phone	Philips PM
Secondary contact Email	Philips PM		
Secondary contact Phone	Philips PM		
Philips Project Manager		Philips Lead installer	Responsibility
Contact name	Philips PM	Contact name	Philips PM
Email address	Philips PM	Phone number	Philips PM
Phone number	Philips PM	Secondary contact name	Philips PM
Secondary contact name	Philips PM	Secondary contact Email	Philips PM
Secondary contact Email	Philips PM	Secondary contact Phone	Philips PM
Secondary contact Phone	Philips PM		
IT department		PRS Lead	Responsibility
Contact name	HCF IT	Contact name	Philips PM
Email address	HCF IT	Email address	Philips PM
Phone number	HCF IT	Phone number	Philips PM
Secondary contact name	HCF IT	Secondary contact name	Philips PM
Secondary contact Email	HCF IT	Secondary contact Email	Philips PM
Secondary contact Phone	HCF IT	Secondary contact Phone	Philips PM

Info & instructions for Philips FSE/Hospital IT (in case of remote connection via ServiceHub)		
ServiceHub	Responsibility	
Static or DHCP	Select Yes or No	HCF IT
IP address	Select Yes or No	HCF IT
Subnet Mask / Subnet prefix length	(should be same for all medical devices)	HCF IT
Default Gateway / Gateway	(should be same for all medical devices)	HCF IT
MAC address		PHILIPS RSE/FSE
Audit trail connection	Select Yes or No	HCF IT
Time Server Connection	Select Yes or No	HCF IT
RSGoC via HSEdge		Responsibility
IP address		PHILIPS PRS lead
Host name		PHILIPS PRS lead
FQDN		PHILIPS PRS lead
Port Number		PHILIPS PRS lead
RSGoC via CMVPN		Responsibility
IP address		PHILIPS PRS lead
Endpoint		PHILIPS PRS lead
RSGoC via Internet		Responsibility
Proxy (optional)	Select Yes or No	PHILIPS PRS lead
IP address (optional)		PHILIPS PRS lead
FQDN (optional)		PHILIPS PRS lead
Port Number (optional)		PHILIPS PRS lead
Username (optional)		PHILIPS PRS lead
Password (optional)		PHILIPS PRS lead
HS Cloud Concentrator (optional)	Select Yes or No	PHILIPS PRS lead
Region		PHILIPS PRS lead
End point		PHILIPS PRS lead
Info & instructions for Bio-Med/Hospital IT (in case of remote connection via PSA/M2M)		
Configure the hospital network for one of the following remote connection types, according to the data supplied on this sheet (open ports in firewall)		
+ VPN connection		
+ VPN connection with source NAT		
+ ISSL connection		
+ ISSL connection with proxy server		


VPN Standard									
Philips Product	Local IP address	Subnet mask	Gateway	Direction	Port	Protocol	Usage	Philips RSN IP range	
Azurion X-ray System	<fill in local IP address>	<fill in local subnet mask >	<fill in local Gateway>	Outbound to	443	HTTPS	For Remote Services Support	https://192.68.49.0/24	
				Inbound from	Ping	ICMP Echo	Reachability check		
				Inbound from	22	SFTP	Azurion R2.1 or higher		
Allura Xper-ray System	<fill in local IP address>	<fill in local subnet mask >	<fill in local Gateway>	Outbound to	443	HTTPS	For Remote Services Support	https://192.68.49.0/24	
				Inbound from	Ping	ICMP Echo	Reachability check		
				Inbound from	21	FTP	Windows XP only		
Interventional workspot	<fill in local IP address>	<fill in local subnet mask >	<fill in local Gateway>	Outbound to	443	HTTPS	For Remote Services Support	https://192.68.49.0/24	
				Inbound from	Ping	ICMP Echo	Reachability check		
Coronary Tools	<fill in local IP address>	<fill in local subnet mask >	<fill in local Gateway>	Outbound to	443	HTTPS	For Remote Services Support	https://192.68.49.0/24	
				Inbound from	Ping	ICMP Echo	Reachability check		
EchoNavigator	<fill in local IP address>	<fill in local subnet mask >	<fill in local Gateway>	Outbound to	443	HTTPS	For Remote Services Support	https://192.68.49.0/24	
				Inbound from	Ping	ICMP Echo	Reachability check		
Philips Hemo	<fill in local IP address>	<fill in local subnet mask >	<fill in local Gateway>	Inbound from	22	SFTP			
VPN with source NAT									
Philips Product		"Source NAT servers - Fixed Philips IP range"				"Source NAT servers - Customer defined symmetrical Source NAT IP"			
Azurion X-ray System		Server Address: https://192.68.49.49 Tunnel Address: https://192.68.49.50 Registration Address: https://192.68.49.140 oneAVS security server: - https://192.68.49.12 - https://192.68.49.29				Server Address : https://62.130.250.49 Tunnel Address : https://62.130.250.50 Registration Address :https://62.130.250.140 oneAVS security server: - https://192.68.49.12 - https://192.68.49.29			
Allura X-ray System									
Interventional workspot									
Coronary Tools									
EchoNavigator									
ServiceHub									
Philips Hemo									
ISSSL standard:									
Philips Product	Local IP address	Direction	Port	Protocol	Philips RSN DNS		Usage		
Azurion X-ray System	<fill in local IP address>	Outbound To	443	HTTPS	https://ws-m2m.prs.healthcare.philips.com https://ta-m2m.prs.healthcare.philips.com https://cam-m2m.prs.healthcare.philips.com		For remote services support		
Allura X-ray System	<fill in local IP address>								
Interventional workspot	<fill in local IP address>								
Coronary Tools	<fill in local IP address>								
EchoNavigator	<fill in local IP address>								
Philips Hemo		No ISSSL Support							
ISSSL with proxy server: Configure "ISSSL standard" and also collect the proxy server settings									
Philips Product	Proxy Server IP address	Proxy Server port number	Proxy Server User name	Proxy Server Password	NTLM (NT Lan Manager)				
Azurion X-ray System	<fill in local Proxy Server IP address>	<fill in local Proxy Server port number >	<fill in local Proxy Server User Name >	<fill in local Proxy Server Password	<NTLM available Yes / No >				
Allura X-ray System									
Interventional workspot									
Coronary Tools									
EchoNavigator									
Philips Hemo									
No ISSSL Support									

(25.0)

X-ray system (Basic)		
Local X-ray system		Responsibility
System Type	E.g. Hospital X	PHILIPS
Serial number		PHILIPS
IP sec enabled	Select Yes or No	HCF IT
HostName	SUITE_PC (always for Azurion)	HCF IT
MAC Address		PHILIPS
IP Address		HCF IT
AE Title	LOCAL_AETITLE	HCF IT
Port Number		HCF IT
AE Title AlluraRIS	AE_ALLURA_RIS	HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Hospital network		Responsibility
Internet protocol used at HCF	IPv4/IPv6	HCF IT
Subnet Mask / Subnet prefix length		HCF IT
Default Gateway / Gateway		HCF IT
Remote service options		Responsibility
Remote Service part of contract	Select Yes or No	PHILIPS PRS lead
RSG available	None/ RSG on-site/ HSEdge/ RSG Cloud/ CMVPN/ RSG Cloud/ Internet	PHILIPS PRS lead
Proactive and Reactive Monitoring	Select Yes or No	PHILIPS PRS lead
Remote Desktop Connection	Select Yes or No	PHILIPS PRS lead
Software Distribution and Installation	Select Yes or No	PHILIPS PRS lead
Coronary Tools		Responsibility
Available in cathlab	Select Yes or No	PHILIPS
Physical location	Control Room/R-cabinet/B-cabinet	PHILIPS
Physical location in (R-/B-)cabinet	Pos1/Pos2/Pos3/Pos4/N.A.	PHILIPS
HostName		HCF IT
MAC Address		PHILIPS
IP Address		HCF IT
AE Title (Storage Node)		HCF IT
AE Title (Storage Commit)		HCF IT
Port Number (Storage Commit)		HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
EchoNavigator		Responsibility
Available in cathlab	Select Yes or No	PHILIPS
Physical location	Control Room/R-cabinet/B-cabinet	PHILIPS
Physical location in (R-/B-)cabinet	Pos1/Pos2/Pos3/Pos4/N.A.	PHILIPS
HostName		HCF IT
MAC Address		PHILIPS
IP Address		HCF IT

Interventional workspot (IW)		Responsibility
Available in cathlab	Select Yes or No	Philips PM
Physical location	Control Room/R-cabinet/B-cabinet	Philips PM
Physical location in (R-/B-)cabinet	Pos1/Pos2/Pos3/Pos4/N.A.	Philips PM
IP sec enabled	Select Yes or No	HCF IT
HostName		HCF IT
MAC Address		PHILIPS
IP Address		HCF IT
AE Title	RSXVimport	HCF IT
Port Number	4110	HCF IT
AE Title IWxRayMod	XVimport	HCF IT
Port Number	3110	HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Network & Services		Responsibility
Windows printer		Responsibility
Physical location	Select Yes or No	HCF IT
Compatible Windows release		HCF IT
IP sec enabled	Select Yes or No	HCF IT
Printer type	e.g. Win10, Win7	HCF IT
HostName		HCF IT
IP Address		HCF IT
Port Number		HCF IT
DICOM printer **		Responsibility
Physical location		HCF IT
Printer type		HCF IT
HostName		HCF IT
IP Address		HCF IT
AE Title		HCF IT
Port Number		HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Structured Report Server 1		Responsibility
Physical location		HCF IT
Model name		HCF IT
HostName		HCF IT
IP Address		HCF IT
AE Title		HCF IT
Port Number		HCF IT
Automatic query schedule date	24 hrs / infinite	HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Certificate name / License		HCF IT


MPPS		Responsibility
Physical location		HCF IT
Model name		HCF IT
HostName		HCF IT
IP Address		HCF IT
AE Title		HCF IT
Port Number		HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Certificate name / License		HCF IT
PPSM IHE compatible	Select Yes or No	HCF IT
Remote system/PACS 1**		Responsibility
Physical location		HCF IT
Model name		HCF IT
HostName		HCF IT
IP Address		HCF IT
AE Title		HCF IT
Port Number		HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Services		
Store AETitle		HCF IT
Store PortNumber		HCF IT
Move AETitle		HCF IT
Move PortNumber		HCF IT
QR AETitle		HCF IT
QR PortNumber		HCF IT
SC AETitle		HCF IT
QR PortNumber		HCF IT
RDSR AETitle		HCF IT
RDSR PortNumber		HCF IT
WorkList Management		Responsibility
Physical location		HCF IT
HostName		HCF IT
IP Address		HCF IT
AE Title		HCF IT
Port Number		HCF IT
Automatic query schedule date	24 hrs / infinite	HCF IT
Secure communication		
Use authentication	Select Yes or No	HCF IT
Use encryption	Select Yes or No	HCF IT
Certificate name / License		HCF IT



1
A1 **SITE LAYOUT (ref.)**
SCALE: 1/4" = 1' - 0"
Required Unistrut Height: 9' - 9 3/8", +1/8" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.

Project Details	Philips Contacts	Project
Drawing Number N-SOU250567_RevC Date Drawn: 8/25/2025 Quote: Q-00552854 Order: 6600762648	Project Manager: Alan Rooks Contact Number: (904) 510-1874 Email: alan.rooks@philips.com Drawn By: Alison Ballargeon	Azurion 7 B20/12 (Cardiac) - Swivel VA Medical Center Gainesville Gainesville, FL Room: C246-1

R-A1.1



2
A1

EQUIPMENT LAYOUT (ref.)

SCALE: 1/4" = 1' - 0"

Required Unistrut Height: 9' - 9 ³/₈" + ¹/₈" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.

R-A1.2

Project Details


Drawing Number
N-SOU250567_RevC
Date Drawn: 8/25/2025
Quote: Q-00552854
Order: 6600762648

Philips Contacts

Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com
Drawn By: Alison Baillargeon

Project

Azurion 7 B20/12 (Cardiac) - Swivel
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1



1
S1

FLOOR & WALL SUPPORT LAYOUT (ref.)


SCALE: 1/4" = 1'-0"

Required Unistrut Height: 9' - 9 3/8", +1/8" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.

NOTE:
ALL SHEETS IN THE "R" SERIES ARE REFERENCES TO VIEWS IN THE MAIN SET AND ARE MEANT TO BE PLOTTED AT 11x17

Project Details	Philips Contacts	Project
Drawing Number N-SOU250567_RevC Date Drawn: 8/25/2025 Quote: Q-00552854 Order: 6600762648	Project Manager: Alan Rooks Contact Number: (904) 510-1874 Email: alan.rooks@philips.com Drawn By: Alison Baillargeon	Azurion 7 B20/12 (Cardiac) - Swivel VA Medical Center Gainesville Gainesville, FL Room: C246-1

R-S1.1



1
S2

CEILING SUPPORT LAYOUT - EQUIPMENT (ref.)

SCALE: 1/4" = 1'-0"

Required Unistrut Height: 9' - 9 ³/₈" , +¹/₈" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.

NOTE:

ALL SHEETS IN THE "R" SERIES ARE REFERENCES TO VIEWS IN THE MAIN SET AND ARE MEANT TO BE PLOTTED AT 11x17

R-S2.1

Project Details

Drawing Number

N-SOU250567_RevC

Date Drawn: 8/25/2025

Quote: Q-00552854

Order: 6600762648

Philips Contacts

Project Manager: Alan Rooks

Contact Number: (904) 510-1874

Email: alan.rooks@philips.com

Drawn By: Alison Baillargeon


Project

Azurion 7 B20/12 (Cardiac) - Swivel

VA Medical Center Gainesville
Gainesville, FL

Room: C246-1

PHILIPS
PROVISIONAL




03.31.2025

Project Details
Drawing Number **N-SOU250567_RevC**
Date Drawn: 8/25/2025
Quote: Q-00552854
Order: 6600762648

R-S2.2

THE INFORMATION IN THIS PACKAGE IS PROVIDED AS A CUSTOMER CONVENIENCE, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS. Philips assumes no liability nor offers any warranty for the fitness or adequacy of the premises or the utilities available at the premises in which the equipment is to be installed, used, or stored.



1 E1 ELECTRICAL PLAN (ref.) --- SCALE: 1/4" = 1'-0" Required Unistrut Height: 9' - 9 ³/₈" , +1" (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut

SCALE: 1/4" = 1'-0"

Required Unistrut Height: 9' - 9 $\frac{3}{8}$ ", + $\frac{1}{8}$ " (2981mm, +3mm) Unistrut height measured from top of Clea floor plate to bottom of Unistrut.

Project
Azurion 7 B20/12 (Cardiac) - Swivel
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com
Drawn By: Alison Baillargeon

Project Details
Drawing Number **N-SOU250567_RevC**
Date Drawn: 8/25/2025
Quote: Q-00552854
Order: 6600762648

R-E1.1

THE INFORMATION IN THIS PACKAGE IS PROVIDED AS A CUSTOMER CONVENIENCE, AND IS NOT TO BE CONSTRUED AS ARCHITECTURAL DRAWINGS OR CONSTRUCTION DOCUMENTS. Phillips assumes no liability nor offers any warranty for the fitness or adequacy of the premises or the utilities available at the premises in which the equipment is to be installed, used, or stored.

PHILIPS
PROVISIONAL

03.31.2025

Point to Point Run Lengths									
Electrician Requirements and Notes									
1. A Greenlee measuring tape (part no. 435) must be provided with raceway (conduit) runs 2. All raceway (conduit) runs must take most direct route point to point. 3. All conduits shall be labeled and identified at each end by run number/destination. 4. For conduit use the minimum bending radius specific to the conduit diameter. All conduit and raceway shall be sweeps.									
	↓	A Raceway (Conduit) supplied/installed by contractor - Philips cables installed by Philips B Raceway (Conduit) supplied/installed by contractor - Philips cables installed by contractor C Raceway (Conduits) and cables supplied and installed by contractor D Raceway (Conduit) existing - cables supplied and installed by Philips E Raceway (Conduit) existing - cables supplied by Philips and installed by contractor F Raceway (Conduit) existing - cables supplied and installed by contractor G Optional equipment, verify with local Philips Service			* P Power (AC) D Power (DC) G Ground S Signal H High Tension C Cooling Hose A Air Supply Hose				
		Raceway (Conduit)		Raceway (Conduit) Quantity	Cable Type (*)	Minimum Raceway (Conduit) Size	Maximum Point to Point Run Length	Special Requirements	(REQUIRED) Electrician to document point to point measurements
		Run No.	From To						
		C 1	ERBThird Party Room Outlets	-	G	-	-	For Injector, Auxiliary Box, Patient Monitoring, Video Networking, etc.	
		C 2	ERB	1	G	3/4"	-	See Sheet "2-ED" for details.	
		A 3	SPME	1	C	2 1/2"	61'	Tube Cooling Hoses.	
		A 4	SPME	1	H	2 1/2"	55'	X-Ray high-voltage cables, the use of 90 deg. ells is not acceptable.	
		A 5	SPMA	1	P/G	1 1/2"	50'	-	
		A 6	SPME	1	P/G	1 1/2"	57'	-	
		A 7	SPMA	1	S	2"	50'	-	
		A 8	SPME	1	S	1"	57'	-	
		A 9	MSAMA	1	P/G	2"	44'	-	
		A 10	MSAMA	1	S	3"	44'	-	
		A 11	MSAMB	1	S	1"	44'	-	
		A 12	PB12ME	1	C	2 1/2"	60'	Tube Cooling Hose.	
		A 13	PB12ME	1	H	2"	62'	X-Ray high-voltage cables, the use of 90 deg. ells is not acceptable.	
		A 14	PB1MA	1	P/G	1 1/2"	59'	-	
		A 15	PB12ME	1	P/G	1 1/2"	63'	-	
		A 16	PB1MA	1	S	3"	59'	-	
		A 17	PB12ME	1	S	1"	63'	-	
		A 18	TVMA	1	P	1 1/2"	60'	-	
		A 19	TVMB	1	P/G	1 1/2"	60'	-	
		A 20	TVMB	1	S	1 1/2"	60'	For FlexVision XL.	
		A 21	TVCY	1	S	1"	60'	-	
		A 22	ATYTV	1	S	1"	60'	-	
		A 23	CYMA	1	S	3"	60'	Conduits to land on wall raceway adjacent to CY. Via raceway.	
		A 24	CYMB	1	S	1 1/2"	60'	Conduits to land on wall raceway adjacent to CY. Via raceway.	
		A 25	CYMA	1	P/G	2"	60'	Conduits to land on wall raceway adjacent to CY. Via raceway.	
		A 26	CYMB	1	S	2"	86'	Conduits to land on wall raceway adjacent to CY. For Philips Workstations. Via raceway.	
		A 27	ATYMA	1	S	2"	47'	-	
		C 28	SN1SN2	1	S	1"	-	For Service Hub. Dedicated conduit and CAT 6 network connection terminated at each end (not allowed with x-ray system cabling).	
		C 29	WLMMA	1	P	PER NEC	PER NEC	X-ray Warning Light Refer to Detail 3-ED	
		C 30	RLMA	1	P	PER NEC	PER NEC	Exam Room Light Refer to Detail 3-ED	

Project	Azurion 7 B20/12 (Cardiac) - Swivel	
	VA Medical Center Gainesville	
Philips Contacts	Gainesville, FL	
	Room: C246-1	
Project Details	Project Manager: Alan Rooks	
	Contact Number: (904) 510-1874	
R-SOU250567_RevC	Email: alan.rooks@philips.com	
	Date Drawn: 8/25/2025	
R-EC1	Quote: Q-00552854	
	Order: 6600762648	

Point to Point Run Lengths									
Electrician Requirements and Notes									
1. A Greenlee measuring tape (part no. 435) must be provided with raceway (conduit) runs									
2. All raceway (conduit) runs must take most direct route point to point.									
3. All conduits shall be labeled and identified at each end by run number/destination.									
4. For conduit use the minimum bending radius specific to the conduit diameter. All conduit and raceway shall be sweeps.									
	A Raceway (Conduit) supplied/installed by contractor - Philips cables installed by Philips							* P Power (AC) D Power (DC) G Ground S Signal H High Tension C Cooling Hose A Air Supply Hose	(REQUIRED) Electrician to document point to point measurements
	B Raceway (Conduit) supplied/installed by contractor - Philips cables installed by contractor								
	C Raceway (Conduits) and cables supplied and installed by contractor								
↓	D Raceway (Conduit) existing - cables supplied and installed by Philips								
	E Raceway (Conduit) existing - cables supplied by Philips and installed by contractor								
	F Raceway (Conduit) existing - cables supplied and installed by contractor								
	G Optional equipment, verify with local Philips Service								
	Run No.	From	To	Raceway (Conduit) Quantity	Cable Type (*)	Minimum Raceway (Conduit) Size	Maximum Point to Point Run Length	Special Requirements	
	G 31	MSA	PHY	1	S	3"	-	For future options (Patient Monitoring). Deleted 12/04/2024.	
	A 32	AUX1	MB	1	S	1 1/2"	86'	Via raceway.	
	A 33	AUX2	MB	1	S	1 1/2"	86'	Via raceway.	
	A 34	AUX3	MB	1	S	1 1/2"	86'	-	
	A 35	AUX4	MB	1	S	1 1/2"	86'	-	
	A 36	AUX5	MB	1	S	1 1/2"	86'	-	
	A 37	AUX6	MB	1	S	1 1/2"	86'	Via raceway.	
	A 38	AUX7	MB	1	S	1 1/2"	86'	-	
	C 39	MSA	TER	1	-	2 1/2"	-	Termination point for for MSA. Deleted 12/04/2024	
	A 40	INT	EPB	1	P	1 1/2"	75'	Conduit opening must be covered if the IntraSight system is planned for future installation.	
	A 41	INT	EPB	1	S	1 1/2"	75'	Conduit opening must be covered if the IntraSight system is planned for future installation.	
	C 42	Power Supply	M3	1	P	3/4"	-	Project electrician to land 120 - 240 VAC, 50 - 60 Hz in M3 terminal block. Light requires: Plus 24, Minus 24, and Ground.	
	A 43	EPB	WR2	1	G	1"	-	Wire to be connected to the ERB.	
	A 44	EPB	CY	1	S	1"	78'	-	
	A 45	EPB	MB	1	S	1"	68'	-	
	C 46	EPB	WR3	1	S	2 1/2"	-	*For 3rd party equipment. Use existing conduit.	
	G 47	EPB	WR3	1	P	2 1/2"	-	*For 3rd party equipment. Use existing conduit.	
	C 48	WR1	WR3	1	P	2 1/2"	-	*For 3rd party equipment to control room. Use existing conduit.	
	G 49	WR1	WR3	1	S	2 1/2"	-	*For 3rd party equipment to control room. Use existing conduit.	
	A 50	IC	MSA	1	S	3"	65'	For Table Mounted Injector.	
	A 51	IC	RIC	1	S	3"	90'	Via raceway.	
	A 52	AVOX1	MB	1	S	1 1/2"	86'	Via raceway.	
	A 53	AVOX2	MB	1	S	1 1/2"	86'	Via raceway.	
	A 54	AVOX3	MB	1	S	1 1/2"	86'	Via raceway.	
	A 55	AVOX4	MB	1	S	1 1/2"	86'	-	
	C 56	Power Panel	CB2	1	P/G	Per N.E.C.	Per N.E.C.	3 AWG minimum phase/ground wire size.	
	C 57	CB2	UPS	1	P/G	Per N.E.C.	Per N.E.C.	3 AWG minimum phase/ground wire size.	
	C 58	UPS	CB	1	P/G	Per N.E.C.	Per N.E.C.	3 AWG minimum phase/ground wire size.	
	C 59	CB	WR2	1	P/G	Per N.E.C.	Per N.E.C.	4 AWG minimum phase/ground wire size. Final connections to terminate at MA terminal block.	
	C 60	SBO	UPS	1	S	Per N.E.C.	Per N.E.C.	For Signaling Box Option. Cable by electrician.	

*
P Power (AC)
D Power (DC)
G Ground
S Signal
H High Tension
C Cooling Hose
A Air Supply Hose

Project

Philips Contacts

Project Details

Azurion 7 B20/12 (Cardiac) - Swivel

N-SOU250567_RevC

R-EC2

VA Medical Center Gainesville

8/25/2025

8/25/2025

Gainesville, FL

Quote: Q-00552854

Order: 6600762648

Room: C246-1

Drawn By: Alison Ballargeon

Point to Point Run Lengths									
Electrician Requirements and Notes									
1. A Greenlee measuring tape (part no. 435) must be provided with raceway (conduit) runs 2. All raceway (conduit) runs must take most direct route point to point. 3. All conduits shall be labeled and identified at each end by run number/destination. 4. For conduit use the minimum bending radius specific to the conduit diameter. All conduit and raceway shall be sweeps.									
	↓	A Raceway (Conduit) supplied/installed by contractor - Philips cables installed by Philips B Raceway (Conduit) supplied/installed by contractor - Philips cables installed by contractor C Raceway (Conduits) and cables supplied and installed by contractor D Raceway (Conduit) existing - cables supplied and installed by Philips E Raceway (Conduit) existing - cables supplied by Philips and installed by contractor F Raceway (Conduit) existing - cables supplied and installed by contractor G Optional equipment, verify with local Philips Service						* P Power (AC) D Power (DC) G Ground S Signal H High Tension C Cooling Hose A Air Supply Hose	
		Raceway (Conduit)			Raceway (Conduit) Quantity	Cable Type (*)	Minimum Raceway (Conduit) Size	Maximum Point to Point Run Length	Special Requirements
		Run No.	From	To					
	B	61	UPS	MA	1	S	Per N.E.C.	Per N.E.C.	UPS communication cable by Socomec.
	C	62	UPS	ST	1	S	Per N.E.C.	Per N.E.C.	
	C	63	CB	ST	1	P	Per N.E.C.	Per N.E.C.	
	G	64	CB2	ST	1	P	Per N.E.C.	Per N.E.C.	
	C	65	TV2	MB	1	S	1 1/2"	86'	For TV2 (Third Party 27" Monitor) Deleted 12/04/2024
	C	66	TV3	MB	1	S	1 1/2"	86'	For TV3 (Third Party 27" Monitor) Deleted 12/04/2024
	A	67	SN1	MB	1	S	1"	26'	
	A	68	MSA	XHA	1	S	2 1/2"	100'	For Xper Hemo App.Use existing 3" conduit.
	A	69	MSA	UPS	1	P	1 1/2"	100'	For Xper Hemo App.Use existing 3" conduit.
	A	70	Boom B	XHA	1	P	1 1/2"	100'	For Xper Hemo App.
	A	71	Boom B	AUX4	1	S	1"	98'	For Xper Hemo App.
	A	72	AVOX2	Boom A1	1	S	1 1/2"	86'	
	A	73	AVOX3	Boom B1	1	S	1 1/2"	86'	

R-EC3

Project Details


Drawing Number
N-SOU250567_RevC
Date Drawn: 8/25/2025
Quote: Q-00552854
Order: 6600762648

Philips Contacts

Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com
Drawn By: Alison Ballargeon

Project

Azurion 7 B20/12 (Cardiac) - Swivel
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1



Legend
Wiring and Connections made by Contractor unless otherwise noted. Isolated from metal casing.
Bond / Safety ground wiring and Connections made by Contractor.
Connections that are connected to metal casing (e.g. outer housing, busbar).
Isolated / Isolated Ground wiring and Connections made by Contractor for equipment grounding conductor. Isolated from metal casing.
Wiring Connection made by Philips
Ground wiring Connection made by Philips
Supply conductors installed together with equipment grounding conductor in the raceway (conduit).

Notes
- Wiring and circuit sizes from source supply must meet Philips regulation requirements and must be determined by contractor to meet building conditions and local codes. 3 Phase + PE (Protective Earth).
- RAP needs to be mounted to the wall by contractor.
- Wire sizes will be dependent on: 1. Lug size of UPS and MA terminal block 2. The tolerance on the allowed nominal mains voltage shall not exceed ±10%. 3. MA CRC Terminal Block can only accept 4 AWG.
- Full UPS assembly, all wiring, installation, Battery Module/Power Module insertion (for modular systems), Seismic brackets (if applicable), and startup in bypass is the responsibility of the electrician unless it is specifically stated otherwise. If electrical work is not completed when commissioning is scheduled there may be a second visit scheduled resulting in delay of timeline and additional cost.
- All Power Cables shall be in separate conduits from control and communication cables.
- ** This equipment may be used for invasive procedures; therefore, the area to be installed is classified as critical care area per NFPA-99 and NFPA-70 (NEC). These documents specify maximum touch voltages and ground impedance in these areas.
- ** The GSSNA specified "Equi-Potential Reference Bar (ERB)" serves as a ground reference for GSSNA equipment. It may also serve as the "Reference Grounding Point" of the room as defined in NFPA 99-3.3.140 for non-Philips Healthcare equipment.

(25.1)

5
ED

ELECTRICAL SCHEMATIC WITH SOCOMEC LLF DETAIL (ref.)

SCALE: N.T.S.

R-ED.5

Project Details
Drawing Number
N-SOU250567_RevC
Date Drawn: 8/25/2025
Quote: Q-00552854
Order: 6600762648

Philips Contacts
Project Manager: Alan Rooks
Contact Number: (904) 510-1874
Email: alan.rooks@philips.com
Drawn By: Alison Baillargeon

Project
Azurion 7 B20/12 (Cardiac) - Swivel
VA Medical Center Gainesville
Gainesville, FL
Room: C246-1

03.31.2025

NOTE:
ALL SHEETS IN THE "R" SERIES ARE REFERENCES TO VIEWS IN THE MAIN SET AND ARE MEANT TO BE PLOTTED AT 11x17